Illumina Genome Analyzer Sequencing Services Microarray Core Facility Huntsman Cancer Institute Room 3363 University of Utah June 2009

Illumina sequencing services performed at the HCI Microarray Core Facility are partially supported by NCI Cancer Center Support Grant P30CA042014 and by institutional support provided through the Huntsman Cancer Institute.

1. Sample Quality Pricing

	Service Category	Description	UofU Pricing	External Pricing
a.	Sample Quality	Bioanalyzer RNA Nano	\$7	\$10.50
b.	Sample Quality	Bioanalyzer RNA Pico	\$8	\$12
c.	Sample Quality	Bioanalyzer Small RNA	\$10	\$15
d.	Sample Quality	Bioanalyzer DNA 1000	\$10	\$15
e.	Sample Quality	Qubit PicoGreen Assay	\$4	\$6
f.	Sample Quality	gDNA gel	\$6	\$9

Notes:

2. Illumina Library Preparation Pricing:

	Service Category	Description	UofU Pricing	External Pricing
a.	Illumina Library Preparation	Genomic DNA	\$385	\$577.50
b.	Illumina Library Preparation	ChIP Seq	\$320	\$480
c.	Illumina Library Preparation	mRNA-seq	\$335	\$502.50
d.	Illumina Library Preparation	Small RNA	\$315	\$472.50

Notes:

¹Genomic DNA provided for Illumina library preps will be qualified using a gDNA gel assay prior to starting the genomic DNA library kit procedure.

¹Recommended concentration of total RNA for Bioanalyzer RNA Nano is 25-500 ng/µl.

²Recommended concentration of total RNA for Bioanalyzer RNA Pico is 50-5000 pg/µl.

³Recommended concentration for Bioanalyzer Small RNA is 50-2000 pg/μl for purified small RNA and 1-100 ng/μl for total RNA.

 $^{^4}$ Recommended concentration of DNA for Bioanalyzer DNA 1000 is 0.1-50 ng/ μ l. This assay covers DNA that ranges from 25-1000 bp. This assay is performed to qualify all Illumina library preparations.

⁵The Qubit PicoGreen Assay is used to determine the concentration ChIP DNA samples.

⁶The gDNA gel assay is used to qualify high molecular weight DNA prior to CGH labeling or Illumina Genomic DNA library preps.

²ChIP DNA provided for Illumina library preps will be qualified using a Qubit PicoGreen assay prior to starting the ChIP-seq library kit procedure.

³RNA provided for Illumina library preps will be qualified using a Bioanalyzer NanoChip gel assay prior to starting the mRNA-seq or the small RNA sample prep kit procedure.

3. Array Capture

	Service/Product Category	Description	UofU Pricing	External Pricing
a.	Array Capture	244K or 1M	\$160	\$240.00
b.	Array Capture	105K or 400K	\$125	
c.	Array Capture	44K or 180K	\$90	
d.	Agilent 244K microarray	all applications/organisms	\$428	\$428
e.	Agilent 105K microarray	all applications/organisms	\$279	\$279
f.	Agilent 44K microarray	all applications/organisms	\$163	\$163
g.	Agilent 1M microarray	all applications/organisms	\$512	\$512
h.	Agilent 400K microarray	all applications/organisms	\$330	\$330
i	Agilent 60K microarray	all applications/organisms	\$194	\$194

Notes:

4. Illumina Genome Analyzer Sequencing Pricing

	Service Category	<u>Description</u>	<u>UofU Pricing</u>	External Pricing
a.	Illumina GA Sequencing	36 Cycle Single Read	\$625	\$937.50
b.	Illumina GA Sequencing	36 Cycle Paired End	\$1155	\$1732.50

Notes:.

¹Genomic DNA for array capture must go through Illumina Genomic DNA library preparation procedure prior to array capture.

Array Capture requires a microarray slide in addition to the array capture service.

³Agilent 244K and 1M microarrays and array capture services using these microarrays must be ordered in multiples of 1.

³Agilent 105K and 400K microarrays and array capture services using these microarray must be ordered in a multiple of 2.

⁴Agilent 44K and 180K microarrays and array capture services using these microarrays must be ordered in a multiple of 4.

¹A 36 Cycle SR sequencing run includes one lane on a flow cell that generates 36 nucleotides of DNA sequence from a single end of each cluster.

²A 36 Cycle PE sequencing run includes one lane on a flow cell that generates 36 nucleotides of DNA sequence from both ends of the DNA clusters within the flowcell lane.